**PX ECOCAST G3**

CE 1250

**PX DENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

BLANC	ALLIAGE DENTAIRE CONVENTIONNEL
WEISS	DENTAL GUSSLEGIERUNGEN
BIANCO	LEGA DENTALE CONVENZIONALE
WHITE	DENTAL CASTING ALLOY
BLANCO	ALEACIÓN DENTAL CONVENCIONAL

% = Au: 3.0, Pd: 23.0, Ag: 65.0, Cu: 6.5, In: 2.0, Ir, Zn: <1%**10.7 g/cm³****990-1050° C****175 HV1****ISO 8891**Fabriqué par / Hergestellt von / Fabricato da /
Made by / Fabricado por**PX DENTAL SA**

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com**PX DENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

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1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX ECOCAST G3 we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX ECOCAST G3 must be cast in a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX ECOCAST G3 casting temperature is 1240° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing : use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a propane/oxygen torch and a pickling paste.

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