**PX GOLD 63****PXDENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

| | |
|----------|--------------------------------|
| JAUNE | ALLIAGE DENTAIRE CONVENTIONNEL |
| GELB | DENTAL GUSSLEGIERUNGEN |
| GIALLO | LEGA DENTALE CONVENZIONALE |
| YELLOW | DENTAL CASTING ALLOY |
| AMARILLO | ALEACIÓN DENTAL CONVENCIONAL |

% = Au: 63.0, Pd: 3.0, Ag: 20.0, Cu: 12.0, Zn: 1.5, Pt, Ir: < 1%

14.3 g/cm³
860-900° C
240 HV1

ISO 1562

Fabriqu e par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds
T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

PXDENTAL SA
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

1. Modelage, alimentation, rev etement

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Instructions for use

1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX GOLD 63 we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX GOLD 63 may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX GOLD 63 casting temperature is 1050° C with a graphite crucible, 1100° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a propane/oxygen torch and a pickling paste.

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