**PX GOLD LF ACTIV**

CE 1250

**PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

JAUNE	ALLIAGE DENTAIRE UNIVERSEL
GELB	DENTAL UNIVERSALLEGIERUNGEN
GIALLO	LEGA DENTALE UNIVERSAL
YELLOW	DENTAL UNIVERSAL ALLOY
AMARILLO	ALEACIÓN DENTAL UNIVERSAL

% = Au: 75.1, Pt: 8.3, Ag: 13.8, Zn: 2.3, In, Ir, Mn,  
Ta: < 1%

16.9 g/cm<sup>3</sup>

980-1070° C

16.0/16.2 µm/m° C

ISO 9693 / 1562

Fabriqué par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

**PX DENTAL SA**

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

**PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

## Mode d'emploi

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX GOLD LF ACTIV we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX GOLD LF ACTIV casting temperature is 1230° C with a graphite crucible, 1270° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min with vacuum at 800° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions. For a better aesthetic use opaque pastes
- PX GOLD LF ACTIV is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

## 8. Soldering

- Pre-bond: use PX SOLD LF 950 solder at 950° C with a flux.
- Post-bond: use PX SOLD 700 solder at 750° C with a flux.

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