**PX GOLD LF****PXDENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

JAUNE	ALLIAGE DENTAIRE UNIVERSEL
GELB	DENTAL UNIVERSALLEGIERUNGEN
GIALLO	LEGA DENTALE UNIVERSAL
YELLOW	DENTAL UNIVERSAL ALLOY
AMARILLO	ALEACIÓN DENTAL UNIVERSAL

% = Au: 75,2, Pt: 9,0, Ag: 9,2, Cu: 3,0, Zn: 2,5, In: 1,0, Ir: <1%**16.6 g/cm³****920-970° C****16.7/17.0 µm / m° C****ISO 9693 / 1562**Fabriqué par / Hergestellt von / Fabricato da /
Made by / Fabricado por**PX DENTAL SA**

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

PXDENTAL SA
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

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1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX GOLD LF we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 700° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX GOLD LF casting temperature is 1120° C with a graphite crucible, 1170° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min with vacuum at 800° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX GOLD LF is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

8. Soldering

- Pre-bond: use PX Sold LF 880 solder at 880° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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