**PX INLAYGOLD PF**

CE 1250

**PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

JAUNE	ALLIAGE DENTAIRE CONVENTIONNEL
GELB	DENTAL GUSSLEGIERUNGEN
GIALLO	LEGA DENTALE CONVENZIONALE
YELLOW	DENTAL CASTING ALLOY
AMARILLO	ALEACIÓN DENTAL CONVENCIONAL

**% = Au: 78.0, Ag: 11.5, Cu: 8.5, Pt, Ir, Zn: <1%****16.0 g/cm<sup>3</sup>****910-950° C****130 HV<sub>1</sub>****ISO 1562**Fabriqué par / Hergestellt von / Fabricato da /  
Made by / Fabricado por**PX DENTAL SA**

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

[www.px dental.com](http://www.px dental.com)**PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

## Mode d'emploi

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX INLAYGOLD PF we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- PX INLAYGOLD PF may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX INLAYGOLD PF casting temperature is 1100° C with a graphite crucible, 1150° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

## 7. Soldering

- Apply the solder PX POSTSOLD 800, yellow, at 800° C using a propane/oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a propane/oxygen torch and a pickling paste.

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