**PX BOND BUDGET****PXDENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACIÓN DENTAL CERAMO-MÉTALICA

%= Au: 15.0, Pd: 52.2, Ag: 20.0, In: 6.0, Sn: 5.5, Ga 1.1,  
Ru : < 1%

12.1 g/cm<sup>2</sup>  
1150-1250° C  
14.1/14.4 µm/m° C

ISO 9693

Fabriqu  par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

**PX DENTAL SA**

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds  
T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

**PXDENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

## Mode d'emploi

### 1. Modelage, alimentation, rev tement

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# Instructions for use

## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX BOND BUDGET we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX BOND BUDGET casting temperature is 1450° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min without vacuum at 950° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX BOND BUDGET is compatible with most known porcelain brands (high fusing).

## 8. Soldering

- Pre-bond: use PX PRESOLD 1110 solder at 1110° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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