**PX ECOCAST PLUS****PXDENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

| | |
|--------|--------------------------------|
| BLANC | ALLIAGE DENTAIRE CONVENTIONNEL |
| WEISS | DENTAL GUSSLEGIERUNGEN |
| BIANCO | LEGA DENTALE CONVENZIONALE |
| WHITE | DENTAL CASTING ALLOY |
| BLANCO | ALEACIÓN DENTAL CONVENCIONAL |

%= Au: 2.0, Pd: 23.9, Ag: 69.0, In: 3.0, Zn: 2.0, Ru : < 1%

10.5 g/cm²
1010-1095° C
180 HV1

ISO 8891

Fabriqu e par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds
T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PXDENTAL SA
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

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Instructions for use

Warning: We advert to our customers that due the high percentage of Silver in this alloy, tarnishing phenomena can occur in mouth.

1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX ECOCAST PLUS we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX ECOCAST PLUS must be cast in a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX ECOCAST PLUS casting temperature is 1295° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing : use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a propane/oxygen torch and a pickling paste.

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