**PX GOLD LF ECO****PX DENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

BLANC	ALLIAGE DENTAIRE UNIVERSEL
WEISS	DENTAL UNIVERSALLEGIERUNGEN
BIANCO	LEGA DENTALE UNIVERSAL
WHITE	DENTAL UNIVERSAL ALLOY
BLANCO	ALEACIÓN DENTAL UNIVERSAL

% = Au: 2.0, Pd: 32.9, Ag: 58.0, Zn: 3.5, In: 1.5,
Sn: 2.0, Ru: < 1%

10.7 g/cm³

1050-1110° C

16.6/16.9 µm/m° C

ISO 9693 / 8891

Fabriqu e par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

PX DENTAL SA
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

1. Modelage, alimentation, rev ttement

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1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX GOLD LF ECO we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- The use of a ceramic crucible is recommended.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX GOLD LF ECO casting temperature is 1310° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10min with vacuum at 800° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX GOLD LF ECO is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

8. Soldering

- Pre-bond: use PX SOLD LF 950 solder at 950° C with a flux.
- Post-bond: use PX SOLD 700 solder at 750° C with a flux.

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