**PX PREMIUM C****PX DENTAL SA**
CH-2300 LA CHAUX-DE-FONDS

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACIÓN DENTAL CERAMO-MÉTALICA

**% = Au: 75.0, Pd: 18.5, Ag: 1.5, In: 2.0, Sn : 2.0, Ir, Ru
Cu, Zn: < 1%**

16.1 g/cm³

1120-1240° C

14.1/14.4 µm/m° C

ISO 9693

**Fabriqué par / Hergestellt von / Fabricato da /
Made by / Fabricado por**

PX DENTAL SA

Rue Numa-Droz 150, 2300 La Chaux-de-Fonds

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.px dental.com

PX DENTAL SA
CH-2300 LA CHAUX-DE-FONDS

Mode d'emploi

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1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX PREMIUM C we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- The use of a ceramic crucible is recommended.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX PREMIUM C casting temperature is 1390° C with a graphite crucible, 1440° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min under vacuum at 950° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM C is compatible with most known porcelain brands (high fusing).

8. Soldering

- Pre-bond: use PX PRESOLD 1100 solder at 1100° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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