**PX PREMIUM D****PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

<b>BLANC</b>	<b>ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE</b>
<b>WEISS</b>	<b>DENTAL AUFBRENNLEGIERUNGEN</b>
<b>BIANCO</b>	<b>LEGA DENTALE PER CERAMICA</b>
<b>WHITE</b>	<b>DENTAL BONDING ALLOY</b>
<b>BLANCO</b>	<b>ALEACIÓN DENTAL CERAMO-MÉTALICA</b>

**% = Au: 80.2, Pt: 4.0, Pd: 13.5, Sn: 2.1, Ir: <1%****17.4 g/cm<sup>3</sup>****1190-1280° C****13.9/14.0 µm/m° C****ISO 9693****Fabriqu e par / Hergestellt von / Fabricato da /  
Made by / Fabricado por****PX DENTAL SA****Rue Numa-Droz 150, 2300 La Chaux-de-Fonds****T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29****www.px dental.com****PX DENTAL SA**  
CH-2300 LA CHAUX-DE-FONDS

## Mode d'emploi

### 1. Modelage, alimentation, rev etement

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX PREMIUM D we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- The use of a ceramic crucible is recommended. Can be casted in graphite crucibles.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX PREMIUM D casting temperature is 1390° C with a graphite crucible, 1440° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min under vacuum at 950° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM D is compatible with most known porcelain brands (high fusing).

## 8. Soldering

- Pre-bond: use PX PRESOLD 1100 solder at 1100° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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