PX EXTRAROCK

PXDENTAL SA

CH.2300 LA CHAUX-DE-FONDS
T. +41 (0)32 922 65 00
F. +41 (0)32 922 65 04

Mode d'emploi

Vos avantages:

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THE IDEAL PRODUCT FOR THE PRODUCTION OF ALL YOUR MODELS.

Working instruction

Your advantages:

- * A product which is both quick and easy to process
- * An excellent fluidity reproducing every detail in the impression
- * A low shrinkage, an excellent dimensional stability and a high degree of accuracy
- * A nonporous material exhibiting a flawless surface finish
- * An exceptional impact resistance, perfect edge stability
- * A choice of colours: the white base can be tinted to choice

Some significant properties:

- * Mixing ratio by weight: 100 A to 20 B
- * Mixing ratio by volume: 100 A to 27 B
- * Onset of hardening: 3-6 min.
- * Shrinkage after 24h: 0.07%
- * De-vest after: 60-90 min.

Compatibility with impression materials

- * Alginates or hydrocolloids: thoroughly eliminate water residue and moisture before pouring the model
- * For the following materials, use release agent PX Extrarock D:
 - Polyether (for example Impregum)
 - Vinyl polysiloxane
 - Thiocacoutchouc elastomer
 - Plastics like burn-out posts
 - All metals including metallic bases of the models etc.

Instructions for use

- 1) carefully clean and dry the impression before use
- 2) Fix the model with a hard catalyst-free silicone or modelling wax to ensure perfect stability.
- 3) In order to avoid surface reactions and bonding when releasing from the mould, apply release agent PX Extrarock D to all surfaces of the impression
- 4) Thoroughly shake component A before opening the tin and pour all the contents into the mixer provided
- 5) Before use, thoroughly mix component A using the mixer provided until perfectly homogeneous
- 6) Pour the desired quantity of component A into the calibrated tumbler using the black calibration marks.
- 7) Colour to choice by adding the dye drop by drop to component A according to the instructions on the bottle (the density of colour is variable and noncritical as the dye is totally miscible)
- 8) Add component B to the component A in the calibrated tumbler up to the corresponding red calibration (the mixing ratio must be observed: i.e. 20% by weight of component B to A if not using the calibration tumbler).
- 9) Thoroughly mix the 2 components using the plastic spatula provided
- 10) Pour the mixture over the impression on a vibrator
- 11) The mould may be placed under a pressure of 2 bars for 15 minutes after pouring
- 12) The onset of hardening begins within 3 to 6 minutes
- 13) It is possible to de-vest after 60-90 minutes

Important points to observe

- * Cleaning and the drying of the impression
- * The application of release agent to the surface of the impression according to the instructions to ensure compatibility with the impression material
- * The complete mixing of component A with the mixer provided
- * Do not overheat the product (use vapour jet for less than 20 sec., do not cure a resin on the PX Extrarock)
- * Before modeling, apply release agent PX Extrarock D on the resulting model.

Precautions related to the product

- * The components A and B are irritating for the eyes, the respiratory tracts and the skin. Handle with care. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. When using PX Extrarock do not eat, drink or smoke.
- * Wear gloves and safety glasses during the handling of the product
- * A safety data sheet is available on request
- * Store the PX Extrarock in a dry room at a temperature of 18 to 20°C
- * Once the containers have been opened the contents should be completely used within 6-9 months
- * Ensure that containers are tightly closed after use
- * Component A: classified under toxicity class 4
- * Component B: classified under toxicity class 3

Products available in this range:

- * Art.N°9101 The basic pack contains: component PX Extrarock A, component PX Extrarock B, calibrated tumblers and spatulas for the mixing of the products
- * Art.N°9102 The mixer: ensures an homogeneous mixture is obtained simply and reliably
- * Art.N°9103 The release agent PX Extrarock D: facilitates the release from the mould and avoids surface reactions
- * Art.N°9104 The dyes: 5 different colours are available