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SOLDER FOR COMBINED WORKS BETWEEN A GOLD BASED ALLOY AND A NI-CR OR CO-CR ALLOY

Composition in %: Au: 80.0, Ni : 12.0, Zn : 8.0

Instructions for use

IMPORTANT DATA

Melting range	840 – 890°C
Working temperature	890°C

Before you start, read the "Important instructions for use and processing of the Dental Alloys". Alloys are products for professional use only.

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

Attention: this product contains nickel which may have allergizing effects or provoke localised reactions. Please consider the use of the product and the particular sensitivity of the patient.

1. The surfaces are held together with wax, resin or composite.
2. The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm².
3. Make the soldering block as small as possible with a sufficient rigidity.
4. Eliminate the wax with boiling water, the resin or composite with an open flame. Degrease the contact surfaces if necessary.
5. Apply a flux to the cold surfaces as well as to the solder.
6. Heat the solder block and the surfaces to be soldered to working temperature (890°C).
7. Apply the solder only when the working temperature is reached.
8. Bench cool the solder and remove the soldering block.
9. Pickle with a usual pickling liquid.

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