PX DENTAL SA

BLANC	ALLIAGE DENTAIRE CONVENTIONNEL
WEISS	DENTAL GUSSLEGIERUNGEN
BIANCO	LEGA DENTALE CONVENZIONALE
WHITE	DENTAL CASTING ALLOY
BLANCO	ALEACION CONVENCIONAL

% = Au: 3.0, Pd:23.0, Ag: 65.0, Cu: 6.5, In 2.0, Ir, Zn: < 1%

10.7 g/cm²
990-1050°C
180 HV1

ISO 22674

PX DENTAL SA

Champs-Montants 16 a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PX DENTAL SA
CH - 2074 MARIN

Mode d'emploi

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Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short or middle span bridges

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX Ecocast G3 may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX Ecocast G3 casting temperature is 1240°C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/ oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a Propane /oxygen torch and a pickling paste.

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