PX DENTAL SA

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACION DENTAL CERAMO-METALICA

% = Au: 2.0, Pd: 78.5, Cu: 6.9, In: 4.5, Sn: 2.0, Ga: 5.5,
Ru, Zn: < 1%

11.4 g/cm³

1150-1270°C

13.9/14.2 µm/m°C

ISO 22674 & 9693

PX DENTAL SA

Champs-Montants 16 a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PX DENTAL SA
CH - 2074 MARIN

Mode d'emploi

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Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures, Cast partial dental frames

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 20-45 min.

3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX EXTRAPALL 2 casting temperature is 1450° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min without vacuum at 950° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX EXTRAPALL 2 is compatible with most known porcelain brands (high fusing).

8. Soldering

- Pre-bond: use PX PRESOLD 1110 solder at 1110° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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