## PX DENTAL SA

|                |                                    |
|----------------|------------------------------------|
| JAUNE CLAIR    | ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE |
| HELLGELB       | DENTAL AUFBRENNLEGIERUNGEN         |
| GIALLO CHIARO  | LEGA DENTALE PER CERAMICA          |
| LIGHT YELLOW   | DENTAL BONDING ALLOY               |
| AMARILLO CLARO | ALEACION DENTAL CERAMO-METALICA    |

% = Au: 84.4; Pt: 7.9; Pd: 4.6, In: 2.6, Ir, Cu: < 1%

18.1 g/cm<sup>3</sup>

1130 - 1190°C

14.2 /14.4 µm/m°C

ISO 22674 & 9693

PX DENTAL SA

Champs-Montants 16 a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PX DENTAL SA  
CH - 2074 MARIN

## Mode d'emploi

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Before you start, read the "Important instructions for use and processing of the Dental Alloys". Alloys are products for professional use only.

## Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

## Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

### 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

### 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

### 3. Crucibles and casting buttons

- The use of a ceramic crucible is recommended.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

### 4. Melting

- For PX PREMIUM S casting temperature is 1340° C with a graphite crucible, 1390° C with a ceramic crucible
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

### 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

### 6. Oxidation

- 10 min under vacuum at 950° C.

### 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM S is compatible with most known porcelain brands (high fusing).

### 8. Soldering

- Pre-bond: use PX PRESOLD 1100 solder at 1100° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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