PX DENTAL SA

JAUNE	ALLIAGE DENTAIRE CONVENTIONNEL
GELB	DENTAL GUSSLEGIERUNGEN
GIALLO	LEGA DENTALE CONVENZIONALE
YELLOW	DENTAL CASTING ALLOY
AMARILLO	ALEACION DENTAL CONVENCIONAL

% = Au: 68.5, Pt:2.3, Pd 4.2 Ag: 12.5, Cu: 10.3, Zn: 2.2
Ir : < 1%

15.0 g/cm²
880-950°C
275 HV1

ISO 22674

PX DENTAL SA

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Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

Indications

Inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures, Cast partial dental frames

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- Qualigold 4 may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For Qualigold 4 casting temperature is 1100° C with a graphite crucible, 1150° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/ oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a Propane /oxygen torch and a pickling paste.

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