PX DENTAL SA

|        |                                    |
|--------|------------------------------------|
| BLANC  | ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE |
| WEISS  | DENTAL AUFBRENNLEGIERUNGEN         |
| BIANCO | LEGA DENTALE PER CERAMICA          |
| WHITE  | DENTAL BONDING ALLOY               |
| BLANCO | ALEACION DENTAL CERAMO-METALICA    |

% = Pd: 59.0, Ag: 28.0, In: 6.5, Sn: 6.0, Ru < 1%

10.7 g/cm<sup>3</sup>

1220-1300°C

14.5 µm/m°C

ISO 22674 & 9693

PX DENTAL SA

Champs-Montants 16 a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PX DENTAL SA  
CH - 2074 MARIN

## Mode d'emploi

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Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

## Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures, Cast partial dental frames

## Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 20-45 min.

## 3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For QUALIPALL 4 casting temperature is 1450° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min without vacuum at 950° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- QUALIPALL 4 is compatible with most known porcelain brands (high fusing).

## 8. Soldering

- Pre-bond: use QUALISOLD P1 solder at 1000-1110° C with a flux.
- Post-bond: use QUALISOLD S2 solder at 735-785° C with a flux.

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