
# Cylinco System Mode d'emploi Gebrauchsanweisung Istruzioni per l'uso

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
# Cylinco System Instructions for Use Modo de empleo

**Produit for Professional use only. Before begin please read "Important Notes on Use and Operation of Dental Attachments"**

## Description

Root post with cylindrical-conical profile. Universal adhesive head, for impression taking and false stump. High precision due to well adapted root canal instruments.

Universally usable posts with adhesive heads for diverging abutment impressions and direct abutments. Flat-headed titanium posts are recommended for abutments and provisional work.

## Indications

Root caps for partial and total dentures, direct and indirect abutments, root post crowns and provisional work.

## Contraindications

When the cavity inside the root canal is too large to achieve continuous contact. Casting of a connector is not possible with non-precious alloys.

### 1. Widening

Prepare the endodontically treated root with the enlarger. The length of the cutting edge of the instruments indicates the maximum depth of the instruments and posts.

### 2. Pre-calibration (Mooser system only)

Work the root canal mechanically with the pre-borer

### 3. Drilling to the final dimension

Calibrate the root canal to the final dimension using the drill and chuck or the contra-angle hand-piece with a very low speed.

### 4. Root canal inlay

An individual root inlay must be prepared to stabilize the abutment or root cap and prevent rotational displacement.

### 5. Adhesive-head posts

Posts made of HMA (palladium-free precious metal alloy) are suitable for universal, repositionable impressions for the fabrication of root caps and direct and indirect abutments.

### 6. Flat milled heads

Titanium posts with flat heads are used for provisional work and composite abutment reconstructions.

### 7. Adhesive bonding

To achieve better bonding between the post and the composite, sandblast the surface of the post with 50 µm aluminum oxide prior to adhesive bonding. (Cementing according to the specific instructions of the manufacturer)

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