
TIMA

Mode d'emploi

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
TIMA

Instructions for Use Modo de empleo

Produit for Professional use only. Before begin please read "Important Notes on Use and Operation of Dental Attachments"

Description

Classic, proven snap-in anchor, low height, with 4 different male parts, for soldering, for connection casting and for casting.

Indications

Removable prostheses on root caps and implants, for example
- total dentures
- intercalated prostheses and mixed intercalated and extended prostheses.

Contraindications

Unilateral free-end dentures without transverse connection

1. Positioning the male part

Using the parallelometer mandrin, position the male part on the flat milled root cap (note the position key) and fix it with adhesive wax. Design the soldering block as small as possible and ensure easy access to the flame.

2. Hardening of the male part

The Au-Pt male parts are self hardening. The male parts welded on the root caps are to be cooled slowly to room temperature.
HMA can be out-casted. After casting, cool slowly to room temperature. When using PMMA burnout resin males, use an alloy with a 0.2% yield strength of at least 500 MPa

3. Polishing

When polishing, be careful not to damage the male part. Protect it with the auxiliary part.

4. Mounting the female part

Before curing, fill the retentive parts with e.g. silicone. The resin must not penetrate between the female and male parts. If there are several anchors, place the female parts parallel to each other. For assembly, position the spacer to allow the resilient solution between the male part between the male part, the root cap and the female part. Remove the spacer disc after the prosthesis has been fabricated. In the mouth, the prosthesis descends under load into the space created by the spacer.

5. Mounting the female part in the mouth

Complete the fabrication of the prosthesis with the auxiliary part and then remove it. The resulting cavity should be slightly larger than the original female part. Drill the dam and place it over the male part, the cap and the gingiva. Place the original female part over the male part. Fill the female part and the retentive parts with e.g. silicone. The resin must not penetrate between the male and female parts. Place a drop of self-curing resin in the cavity and insert the prosthesis. After the resin has cured, remove the prosthesis. If necessary, fill the remaining cavity with fresh resin and cure in a pressure pot if possible. Remove the excess resin first.

6. Temporary solution with the dummy

The auxiliary part can replace the original female part and serve as a temporary solution.
If there are several anchors, the force required to remove the prosthesis can be minimized. Before mounting in the resin, lightly grind the outer surfaces of the auxiliary part.

7. Replacing the spring ring

Insert the screwdriver into the closing screw, unscrew it by turning to the left and remove it. Replace the spring and tighten the locking screw with the screwdriver. Simplify the assembly of the closure screw and the spring in the prosthesis by screwing them into an empty female part (the spring is held by the screwdriver)

8. Activation

Insert the screwdriver into the locking screw, unscrew it by turning to the left and remove it. Tighten the locking screw with integrated hard spring

9. Removal of the female part

Remove the locking screw and insert the puller into the empty female part. Heat the extractor, hold it with forceps and pull the socket out of the denture resin.

10. Relining

Root caps: To make the model, place the transfer pin in the impression or in the female part of the prosthesis, then make the model.
Implant work: Instead of using a Timas transfer pin, use a suitable replica or the implant analog

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